

THE SPECIALIST FOR MICROSURGERY

Create your perfect sitting position - with the three-dimensionally adjustable armrests and the patented seat upholstery. Keep your back healthy. Remain sterile and simply adjust your height via the footswitch. High-quality components make the BALANCE ADVANCE a long-lasting working tool.



MODELS

BALANCE ADVANCE B*

ART. NO. H2.12.12
with height adjustment
via foot activation
& a brake

BALANCE ADVANCE

ART. NO. H2.11.12
with foot-activated
height adjustment

*BRAKE

OR CHAIR

BALANCE ADVANCE


FEATURES

- » Three-dimensionally adjustable armrests
- » Stable and backlash-free fixation of armrest position
- » Height adjustment from 590-790 mm
- » Central brake via foot activation
- » Adjustable backrest with lordosis support
- » Patented seat shape with adjustable seat angle
- » Comfort undercarriage with antistatic wheels
- » Quality manufactured with high-quality materials
- » Attractive design

TECHNICAL DATA

5 ball-bearing polyurethane castors, one of them conductive	Ø 65 mm
Height adjustment	590 - 790 mm
Inclination of backrest	+6° / -16°
Tilting range of backrest cushion	25°
Height adjustment of backrest	100 mm
Inclination of seat	+5° / -8°
Seat dimensions (W x H x D)	460 x 60 x 420 mm
Weight	20.3 kg
Maximum load	120 kg

| THREE-DIMENSIONAL ADJUSTABLE



| HIGH-QUALITY MATERIALS



| ADJUSTABLE BACKREST

