

MODELS

BALANCE B*

ART. NO. H1.10.11
with height adjustment
via foot activation
& a brake



BALANCE

ART. NO. H1.00.11
with foot-activated
height adjustment

MODELS

BALANCE PLUS B*

ART. NO. H1.11.11
with height adjustment
via foot activation
& a brake



BALANCE PLUS

ART. NO. H1.01.11
with foot-activated
height adjustment



*BRAKE

OR CHAIR

BALANCE

FEATURES

- » Flexible height adjustment from 480 to 680 mm
- » Patented seat shape with comfortable adjustable seat angle
- » Adjustable backrest with lordosis support
- » High-quality undercarriage with braking, antistatic wheels and hygienically enclosed, enlarged castors
- » Easy-to-operate foot activators for adjusting the seat height without loss of sterility

Optional

- » Column elevation 100 mm (Article no. E.001558)

OR CHAIR

BALANCE PLUS

| NUMEROUS SEAT ADJUSTMENTS



| BACK-FRIENDLY



| BRAKING FUNCTION (OPTIONAL)

