

MODELS

COMFORT F/B*

ART. NO. H1.10.00
with height adjustment
via foot activation and brake

COMFORT H

ART. NO. H1.00.01
with manual
height adjustment

COMFORT SWING F/B*

ART. NO. H1.12.00
with height adjustment
via foot activation and brake

COMFORT SWING H

ART. NO. H1.02.01
with manual
height adjustment

COMFORT SWING F

ART. NO. H1.02.00
with foot-activated
height adjustment

COMFORT F

ART. NO. H1.00.00
with foot-activated
height adjustment

*BRAKE

THE ALL-ROUNDER

The simply designed and easy-to-operate work stool provides valuable service in practices and ORs. An inexpensive and sturdy seating option.



OR STOOL

COMFORT

FEATURES

- » Comfortable and high-frequency-welded seat padding for reliable hygiene and easy cleaning
- » Undercarriage with large, enclosed, antistatic and therefore also hygienic castors
- » Flexible height adjustment from 540 to 725 mm (Comfort H) or from 510 to 710 mm (Comfort F, F/B)

Optional

- » Column elevation 100 mm (Article no. E.001558)


OR STOOL

COMFORT SWING

| MANUAL ADJUSTMENT



| FOOT ACTIVATION



| COMFORT SWING SYSTEM

