

NEW
THE PRO FOR
OPHTHALMOLOGY

ART. NO. Y1.02.12

OPERATING TABLE

GENIUS EYE PRO

FEATURES

- » Telescopic head section with horseshoe-shaped headrest
- » Memory function with 8 programs
- » Battery pack operation 24V DC (for up to 2 weeks)
- » Remote control
- » Telescopic column 595-885 mm
- » Stainless steel construction
- » Push handles at head and foot end
- » 250 kg patient load
- » Additional mounting rails
- » Extra thick padding

- Optional:
- » Dual telescopic column 535-965 mm
 - » Attachment of various armrests
 - » 230 V mains operation
 - » Footswitch

UNDERCARRIAGES



4-EASY-DRIVE - GENIUS


- » 4 guide castors, with connecting axle & lift-brake system
- » Intelligent directional guide
- » Maximum legroom
- » Smooth-running wheels 120 mm



5-WHEEL DRIVE

- » Selective 5th wheel
- » Best riding comfort for frequent users
- » Maximum legroom
- » Smooth-running wheels 120 mm

| FLEXIBLY ADJUSTABLE ARMRESTS



| PUSH HANDLES AT HEAD END



| ADDITIONAL MOUNTING RAILS

