



## THE EXPERT FOR OPHTHALMOLOGY

The GENIUS EYE was optimized perfectly for ophthalmology and the associated surgical procedures.

ART. NO. Y1.01.12

### OPERATING TABLE

# GENIUS EYE

#### FEATURES

- » Telescopic head section with horseshoe-shaped headrest
- » Memory function with 8 programs
- » Battery pack operation 24V DC
- » Remote control
- » Telescopic column 580-870 mm
- » Integrated armrests
- » Stainless steel construction
- » Push handles at head and foot end

#### Optional:

- » Dual telescopic column 520-950 mm
- » 230 V mains operation
- » Footswitch

#### UNDERCARRIAGES



#### 4-EASY-DRIVE - GENIUS

- » 4 guide castors, with connecting axle & lift-brake system
- » Intelligent directional guide
- » Maximum legroom
- » Smooth-running wheels 120 mm



#### 5-WHEEL DRIVE

- » Selective 5th wheel
- » Best riding comfort for frequent users
- » Maximum legroom
- » Smooth-running wheels 120 mm

| INTRAOPERATIVELY ADJUSTABLE HEADREST

| PUSH HANDLES AT HEAD END

| INTEGRATED ARMRESTS

