

A GENIUS IN ITS SEGMENT

If you need a sturdy base with premium quality, then GENIUS is top choice. The economical solution from BRUMABA.



ART. NO. G1.01.12

OPERATING TABLE

GENIUS

FEATURES

- » Telescopic column GENIUS from 580-870 mm
- » Telescopic head section with horseshoe-shaped headrest
- » Memory function with 8 programs
- » 230 V mains operation
- » Back section extension for length adjustment
- » Remote control
- » Stainless steel construction

Optional:

- » Battery pack operation 24V DC
- » Dual telescopic column GENIUS from 520-950 mm
- » Footswitch

UNDERCARRIAGES



2-GUIDE-DRIVE "CL" - GENIUS

- » Two directionally guided castors with levered undercarriage (head end)
- » Two guide castors, with connecting axle and lift-brake system

4-EASY-DRIVE - GENIUS

- » 4 guide castors, with connecting axle & lift-brake system
- » Intelligent directional guide
- » Maximum legroom
- » Smooth-running wheels 120 mm

| EXCELLENT MOBILITY



| PERFECT ACCESS



| ADJUSTS TO YOUR POSITION

