

## BEST IN CLASS

The special lateral tilting feature by 17° and extreme height adjustment allow back-friendly working during sitting or standing.



ART. NO. 02.01.31

### OPERATING TABLE

# PRIMUS

#### FEATURES

- » Telescopic column 520-950 mm
- » Lateral tilting
- » Memory function with 8 programs
- » Telescopic head section with horseshoe-shaped headrest
- » Lever and steering castor undercarriage
- » 24V DC battery pack operation
- » Remote control
- » Back section extension for length adjustment
- » Stainless steel construction

#### UNDERCARRIAGES



#### 2-GUIDE-DRIVE - PRIMUS

- » Two directionally guided castors with levered undercarriage (head end)
- » Two guide castors, with connecting axle and lift-brake system



#### 4-EASY-DRIVE - PRIMUS

- » 4 guide castors, with connecting axle & lift-brake system
- » Intelligent directional guide
- » Maximum legroom
- » Smooth-running wheels 120 mm

#### | COMFORTABLE ENTRY HEIGHT



#### | LATERAL TILTING



#### | BEST POSSIBLE ACCESS

