



# THE UNIVERSAL ALL-ROUNDER

The name says it all: the VARIUS covers virtually all surgical disciplines.

ART. NO. 02.02.31

## OPERATING TABLE

# VARIUS

### FEATURES

- » Dual telescopic column 610-1,040 mm
- » Lateral tilting
- » Memory function with 8 programs
- » Divided foot section
- » Telescopic head section with horseshoe-shaped headrest
- » Back section extension for length adjustment
- » Lever and steering castor undercarriage
- » 24V DC battery pack operation
- » Remote control

### UNDERCARRIAGES



#### 2-GUIDE-DRIVE "OR" - VARIUS

- » Two directionally guided castors with levered undercarriage (head end)
- » Two guide castors, with connecting axle and lift-brake system



#### 4-EASY-DRIVE - VARIUS

- » 4 guide castors, with connecting axle & lift-brake system
- » Intelligent directional guide
- » Maximum legroom
- » Smooth-running wheels 120 mm

### | BEACH CHAIR POSITION



### | PERFECT ACCESS



### | X-RAY CAPABILITY

